

OpeMed Bea

innovative solutions



BEA

A versatile, simple to use product that can be used for showering, toileting, dressing or simple and quick transfers

OpeMed launches the BEA in UK following success in Europe

The BEA services a whole range of manual handling needs from showering, toileting, dressing whilst minimising the number of transfers.

The BEA's fully stainless steel design means it is the only product that can be used for showering, and the shape allows use with wider patients.

Its simplistic design aids training, and halves the number of transfers from a rotastand, as the patient can be transported on the BEA.

Options include a commode seat, extra wide model, reverse leg model (for specialist powered chairs) and an adjustable base



Keeping the "nose above the toes", the patient is encouraged to stand in an ergonomic way for dressing or washing



The BEA ensures the patient is perched in a safe upright position, making standing again quick and effortless



The seat is quickly positioned on the retractable supports and the back rail is reversed to provide security when transferring or lowering onto a chair or WC



The legs of the BEA fit around a WC and between a wheelchair or chair. The seat is at a height to fit over the bed. The double castors prevent sinkage and allow optimum manoeuvrability

OpeMed (Europe) Ltd

t: +44 (0)1252 758858

f: +44(0)870 161 9681

e: sales@opemed.net

19-20 Riverside Park, Farnham, Surrey. GU9 7UG. United Kingdom

www.opemed.net

