


**DIMENSIONS FOR TILT-IN-SPACE SHOWER CHAIR AND CRADLE**

## Technical Specifications for Standard Product

### Adult T.I.S. Cradle

Seat Width	Overall Width	Seat Height	Overall Height	Seat Depth	Overall Depth	Overall Depth When Tilted	SWL	Weight
								
480mm	650mm	AJD	1320mm	430mm	1100mm	1450mm	150kg	17kg

### Paediatric T.I.S. Cradle

Seat Width	Overall Width	Seat Height	Overall Height	Seat Depth	Overall Depth	Overall Depth When Tilted	SWL	Weight
								
450mm	610mm	AJD	1190mm	380mm	990mm	1190mm	150kg	15kg

### Adult T.I.S. Chair

Seat Width	Overall Width	Seat Height	Overall Height	Seat Depth	Overall Depth	Overall Depth When Tilted	SWL	Weight
								
480mm	650mm	AJD	1370mm	450mm	1100mm	1450mm	150kg	20kg

### Paediatric T.I.S. Chair

Seat Width	Overall Width	Seat Height	Overall Height	Seat Depth	Overall Depth	Overall Depth When Tilted	SWL	Weight
								
450mm	610mm	AJD	1190mm	380mm	990mm	1190mm	150kg	18kg

Special Sizes available

### OPTIONS

Neoprene Harness, Commode Pan, Padded Seat, Detachable Wedges