



Solutions for Transfer and Care

## Changing Places Arrives at the Prestigious Guildhall, Swansea



The Guildhall has undergone a huge amount of regeneration work to modernise the impressive multipurpose entertainment venue whilst keeping and restoring many of the original features.

Incorporating accessible features has been very important for this project and the venue boasts a number of wheelchair friendly entrances and exits that are all clearly marked. There is also the added benefit of a dropped bar area for wheelchair users.



# aveso

official sponsor to changing places campaign  
[www.aveso.co.uk](http://www.aveso.co.uk)

## Case Study

Changing Places Facility

The Guildhall at Brangwyn Hall,  
Swansea, Wales



OpeMed were called upon to help the Guildhall with their Changing Places project and guide them through each stage including equipment selection

An initial site visit took place by our Product Specialist including CAD drawings of the room layout and a quotation for the required equipment.

Once the preliminary work and electrics were in place, a qualified OpeMed engineer visited the site to ensure everything was correct and ready for installation in accordance with our pre installation requirements.

The specialist equipment chosen for the Guildhall Changing Places toilet included the OpeMed bariatric Changing Table and the OT200 Compact Ceiling Hoist on an XY to XY system for full room coverage.

*OpeMed are the official co-sponsor of the Changing Places Campaign as Aveso, visit [www.aveso.co.uk](http://www.aveso.co.uk) to view our installs.*



Solutions for Transfer and Care

## Case Study

Changing Places Facility  
The Guildhall at Brangwyn Hall,  
Swansea, Wales

### The Equipment Installed

The OpeMed Bariatric Changing Table is fully height adjustable at the touch of a handset and folds away neatly against the wall when not in use.

The Changing Table stretcher is made of a special high tension tear resistant material with drain holes and integrated water tray for ease of showering.

The OpeMed OT200 Compact Ceiling Hoist chosen for The Guildhall, Swansea was a 2 Function Ceiling Hoist chosen for its full room coverage and ease of use from the simple handset control.

The hoist allows users to be transferred from the toilet to Changing Table or wheelchair with ease and dignity.

Changing Places are life-changing facilities for many thousands of people who require extra care and are unable to use standard accessible toilets comfortably.



The Guildhall in Swansea opened in 1934 and plays a central role within the lives of local people, now including more disabled visitors



Without Changing Places, people with profound and multiple learning disabilities as well as other serious impairments such as spinal injuries & muscular dystrophy often need additional support from a ceiling hoist or changing bench.

*For more information about campaigning or installing a Changing Places please contact us on [info@aveso.co.uk](mailto:info@aveso.co.uk)*

Hoists

Ceiling Hoists

Changing Tables

Showering Tables

Shower Chairs

Bathing

OpeMed (Europe) Ltd  
Clockhouse, Farnham  
Surrey. GU9 7UD  
UNITED KINGDOM

[www.opemed.net](http://www.opemed.net)

T: +44 (0)1252 758858  
E: [sales@opemed.net](mailto:sales@opemed.net)  
F: +44 (0)870 1619681