


OpeMed Install Assisted Living Equipment into Oxfordshire Home

In a recent Case Management adaptation project, OpeMed led the renovation to transform a house into a well adapted home to meet the needs of the disabled user and their family.


Working with Case Managers, OpeMed supplied a Through Floor Lift, full XY Ceiling Tracking upstairs and downstairs, Assisted Bath, Shower Trolley, Assisted WC and Height Adjustable Wash Basin to adapt both the upstairs and downstairs of the house.



➤ **Downstairs Assisted Bathroom:** Hi-lo Wash Basin, assisted W/C and hydraulic Shower Trolley



➤ **Downstairs Assisted Bathroom:** Height adjustable Wash Basin. Through Floor Lift from ground to first floor




➤ **Upstairs Assisted Bathroom:** OT200 Ceiling Hoist, Bath with Changing Table, Assisted W/C and Wash Basin



➤ **Therapy Room:** OT200 compact Ceiling Hoist on XY to XY frame for full-room coverage

Installed Equipment



The OT200 Compact Ceiling Hoist

Ideal for comfortable and easy transfer from Shower Trolley to W/C or Changing Table to Hi-Lo Bath in this case. The Ceiling Hoist is ideal for daily use as it works on a gravity powered system which charges during use. The 6 hook spreader bar is designed for easy sling adjustment and will pick a user up from as low as the floor.

The OpeSpray Shower Trolley

The height adjustable hydraulic Shower Trolley is suitable for both showering, transferring and changing. The Shower Trolley height is adjusted via foot pedal and 4 degree tilt to the head or foot is operated via lever. The Shower Trolley can be transferred between

rooms or left stationary depending upon the user's needs.

Bath with Integrated Changing Table -

The OpeMed Hi-Lo Bath with Changing Table offers a space saving option for transferring and drying users all in one space. When not in use, the Changing Table folds in half to one side allowing the bath to be used. The two in one bath works well with the OT200 Ceiling Hoist as a time saving combination of equipment.

Assistive W/C & Wash Basin –

The OpeMed assistive W/C and Wash Basin options vary depending upon the needs of the user. Our equipment adjusts to the correct height of the user, can be supplied in unique shapes and colours as well as storing personalized preferences for hygienic washing and drying.

