



## Rio Height Adjustable Bath

The Rio bath (with optional patented Neatfold stretcher as shown) is made using the highest quality materials including the twin actuators, which make this one of the quietest height adjustable baths available today.


This beautifully designed and finished bath will enhance any bathroom environment, allowing the client the dignity of taking a bath in the comfort of their own home, while also allowing the carer to have the pleasure and convenience of bathing the client at a safe and comfortable working height, helping to reduce the risk of back strain. As with all OpeMed disabled baths, the Rio will comfortably accommodate all types of internal body supports.

### Unique Features

- Optional Universal Footboard
- Height adjustable
- Smooth and quiet operation
- Handset control
- Additional Neatfold Changing Table
- Available with Music, Lights and Air Spa options



Air Jets, Lights and Music Therapy System




Additional Option of Neatfold Changing Table.




Sturdy base for additional stability

**Product Specifications** Rio Height Adj. Bath

# Specifications

Model	Rio Height Adjustable Hydrotherapy Bath
Function	Height Adjustable Bathing, Hydrotherapy Bathing, Assisted Bathing
Dimensions	1695 x 755mm
Height Range	700mm -1000mm Height Range,
Weight Capacity	250kg safe working load
Options	Range of Taps including Thermostatic and shower, Neatfold Changing Stretcher System, Bluetooth Music System, LED Lights, Warm Air Spa, Universal Footboard, 3 Stage backrest
Application	Nursing Home, Hospital, Domestic, Hospice
LED Light Colours	

# Dimensions



Rio/Kiva Bath with Neatfold Stretcher

