






# Spreader Bar Options



## Standard Spreader Bars

All standard spreader bars have a **200kg or 250kg** safe working load.

	Type	Sling Type	Features
	2 Point Spreader Bar	Loop	Included as standard with the OT200 and OT200R models, the 2 Point Spreader Bar is an easy to use system, allowing users and carers to seamlessly hook looped slings onto either of the points. Spreader bar also comes with padded central bar for additional protection. Code: OT200-32000331
	Manual Cradle	Clip	The Manual Cradle is a self balancing cradle that moves with the body. As the user is hoisted from a reclined position to a seated position, the Manual Cradle naturally follows the user's movement helping to reduce the risk of skin shear to the user.  The Manual Cradle has been designed to provide the carer with a full view of the user offering the carer the chance to remain engaged with the user during the transfer. Code: 32000020
	Powered Cradle	Clip	The Powered Cradle is an assisted cradle powered by an integrated actuator. As the user is hoisted from a reclined position to a seated position, the Powered Cradle moves with the body and the change of angle, helping to reduce the risk of skin shear to the user. Code: 32000140
	X-Frame 4 Point Spreader Bar	Loop	The X-Frame 4 Point Spreader Bar is designed to spread the users weight across four individual support arm for each corner of the body. Having this feature enables the user to be easily supported without the sling becoming bunched above the users head. Code: 32000630
	Quick Release System	Clip / Loop	Compatible with all standard spreader bars, the OT200 Quick Release System allows a no-tools quick change from one spreader bar to another. Code: 3210550

## Bariatric Spreader Bars

All bariatric spreader bars have a **400kg** safe working load.

	Type	Sling Type	Features
	Powered Bariatric Cradle	Clip	The Powered Bariatric Cradle is an assisted cradle powered by an integrated actuator designed to work with up to 400kg. As the user is hoisted from a reclined position to a seated position, the Powered Cradle moves with the body supporting the user and reducing the risk of skin shear. Code: 32000200
	X-Frame Bariatric 4 Point Spreader Bar	Loop	The X-Frame Bariatric 4 Point Spreader Bar is designed support up to 400kg of weight spread across four individual support arms for each corner of the body. This enables the user to be supported without the sling becoming bunched above the users head. Code: 32000630