



Comfortline Height Adjustable Bath

The OpeMed Comfortline is a freestanding, aesthetically pleasing, height adjustable bath. Carefully designed to ensure a relaxing and therapeutic bathing experience, the comfortline at a practical working height for the carer. The bath can be positioned in almost any central bathing area or bathroom and when used in conjunction with a hoist, the Comfortline is easily accessible from either side or either end of the bath.


The Comfortline can come with additional hydrotherapy options for the user. The bath includes key features such as rotating sensory lights, air spa and a music sound wave system to create the desired bathing environment.

Unique Features

- Range of Taps incl. TMV Temperature Controls
- Options of Air Spa, Lighting and Music
- Easy operation 6 hook spreader bar for easy sling adjustment
- Detached – No wall or floor fixings
- Available in Choice of RAL Bespoke Colour
- Ergonomic and Attractive Design Tub Design
- Large Internal Bathing Space



Adjustable Foot Plate




Customisable Hydrotherapy Experience



Available with Lights, Air Jets and Music



Available in bespoke colour




Range of taps available



Ergonomic Tub Design

Comfortline Bath
Product Specifications

Specifications

Model	Comfortline Height Adjustable Hydrotherapy Bath
Function	Height Adjustable Bathing, Hydrotherapy Bathing, Assisted Bathing
Dimensions	2080mm x 780mm
Weight Capacity	150kg safe working load
Options	Autofill with 2 pre-set levels, Optional RAL colours, Fixed height bath (Ecoline), Air-spa system with 12 integrated jets, Mp3 music sound wave system with 8 intergrated speakers, LED colour light system with 6 options, Range of Taps incl. TMV Temperature Controls, Available in Choice of RAL Colour, Removable Changing Stretcher.
Application	Nursing Home, Hospital, Domestic, Hospice
LED Light Colours	

Dimensions

