

OpeMed Track Span Details



90mm Track Span

- Up to **2.6m** when wall mounted by wall post, wall bracket or gantry.
- When using a ceiling bracket, track is supported from above and can be used in multiples through use of a joining bracket.



140mm Track Span

- Up to **4.5m** when wall mounted by wall post, wall bracket or gantry.



180mm Track Span

- Up to **6.0m** when wall mounted by wall post, wall bracket or gantry.

Please Note: Spans mentioned above are for Ceiling Hoists rated 200kg safe working load. If the Hoist is to be down rated to 120kg, the span may be increased but the office must be consulted prior to any agreement with clients for a technical consultation.

For further information on our other tracking systems, please call
+44 (0)1252 758 858

or visit
www.OpeMed.net/product/ceiling-track/

