



OpeTrack Gantry System

The OpeMed Gantry is a portable track system ideal in situations where permanent ceiling track isn't necessary. Compatible with fixed and portable OpeMed Hoists, OpeTrack offers a temporary solution or an alternative hoisting system where installing track into the ceiling is not possible.

Ideal for use with the Monarch Portable Hoist, the OpeTrack Gantry has a weight capacity of up to 200kg. Unobtrusive and space saving against other Gantry systems; it has one of the thinnest base sizes available on the market which greatly reduces tripping hazard for those working around them and with an easy-to-install system, can be tailor made to the room.

Unique Features


- Weight Capacity of 200kg
- Extra low profile base to reduce trip hazard
- Easy assembly, no drilling required
- Easy and fast to tailor on-site



Thin base to reduce trip hazard



Easy to assemble with a single Allen Key



Suitable for fixed & portable OpeMed Ceiling Hoists



Height and length can be easily tailored on-site



Available in subtle off-white colour



Ideal for the Monarch Portable Ceiling Hoist

OpeTrack Gantry
Product Specifications

Specifications

Model	Ope-Track Gantry System
Function	Patient Lifting, Room-to-Room Transfer, Ideal for temporary installations
Weight Capacity	200kg
Dimensions	X + 240mm x 2438mm
Application	Holiday homes, Hotels, Hospice, Rental Property, Domestic, Nursing Home, Hospital

Dimensions

