

Case study - Australia:

Hospital Becomes State-of-the-Art with Hoist Systems



Hospital Becomes State-of-the-Art with Multi-room Hoist Systems

The Hospital

Working closely with our Australian partners, the state-of-the-art Bendigo Hospital room layouts have been equipped with Straight and Curved Track, XY (H) Full Room Covering System incorporating both fixed (OT200) and Portable hoists system throughout Bendigo Hospital dedicated healthcare facilities.

The areas implemented are MRI (Magnetic Resonance Imaging) Suite, Emergency, Holding/Recovery, Ultrasound, Short Stay, Aboriginal Services, Integrated Cancer Centre, Cardiac Recovery, Plaster, Birthing Suite, Outpatients, Women Health, Children & Adolescent, Rehab, Surgical, Medical Rooms and Mortuary.

Working with our dealership in Australia, we managed and accounted the entire project from the end-to-end process which includes designing, installation, commissioning and training.

- 3 therapy rooms with Fixed Hoist Systems
- 38 multi-room bedroom to bathroom Fixed Hoist Systems
- 21 rooms with Fixed Hoist Systems mounted onto Straight Tracking
- 110 rooms with Portable Hoist Systems mounted onto Straight Tracking

OpeMed works closely with Architects, Developers, Builders through global dealers to ensure project objectives meet the countries national standards.



The Installation

Straight Track System: The single rail is positioned perpendicular to the bed, providing coverage from bed to chair. This rail system is suitable to connect on a portable hoist as and when required by a patient.



XY Gantry Ceiling Rail System: The fixed hoist system allows for complete mobility from a bed to the bathroom and vice versa.



Curved Track System: The Curved Track allows a smooth glide when transferring a patient from bed to the bathroom and providing assistance for patient handling.

XY (H) Full Room Covering System With Access To Bathroom: The XY Gantry Ceiling Rail System is connected using a Gate Assembly, allowing caregiver to move a client directly from Bed/Chair to Bathroom.

XY (H) Full Room Covering System in Birthing Suite: The versatility of the ceiling hoist system provides patient with mobility assistant after performing surgery.

XY (H) Full Room Covering System in Gym Therapist room: The fixed hoist mounted on an XY Gantry System assist patients with mobility and speedy recovery.

XY (H) Full Room Covering System in Smaller Suite: The XY Gantry System can be fitted in a smaller suite, covering a full-length bed on XY distance, transfer from bed to a chair and pick up from the floor.




International Projects

OpeMed work closely on an international scale with dealerships placed throughout Europe, Asia, the Middle East and Australasia to provide OpeMed based solutions for hospitals, long-term care, rehab centres and home care.

Alongside our diverse range of solutions we offer a range of tools and support for our international colleagues, visit our website to learn more.

www.opemed.net/our-customers/international-clients-patient-transfer-equipment/



For more information, contact the **OpeMed Team** and call **+44 (0) 1252 758858**

www.opemed.net